Albert-Ludwigs-Universität Freiburg, Institut für Informatik

Prof. Wolfram Burgard, Alexander Schaefer

Lecture: Robot Mapping Winter term 2016/2017

General Notice

Sheet 1

Topic: Octave
Due: November 24, 2016

- The exercises should be solved in groups of two students.
- We will be using Octave for the programming exercises. Octave is a command line program for solving numerical computations. The software is an open-source imitation of its commercial counterpart MATLAB and mostly compatible with it. The software is freely available from www.octave.org. It is compatible with Linux, Mac OS, and Windows. Install Octave on your system in order to solve the programming assignments. A quick guide to Octave is given in the Octave cheat sheet which is available on the website of this lecture.

If you rather want to use MATLAB, which is freely available to all students of the University of Freiburg, you are responsible for adapting the example code, if necessary.

Exercise 1: Get familiar with Octave

The purpose of this exercise is to familiarize yourself with Octave and learn basic commands and operations that you will need when solving the programming exercises throughout this course.

Go through the provided Octave cheat sheet and try out the different commands. As pointed out in the sheet, a very useful Octave command is help <function>. Use it to get information about the correct way to call any Octave function.

Exercise 2: Implement an odometry model

Implement an Octave function to compute the pose of a robot based on given odometry commands and its previous pose. Do not consider the motion noise here.

For this exercise, we provide you with a small Octave framework for reading log files and to visualize results. To use it, call the main.m script. This starts the main loop that computes the pose of the robot at each time step and plots it in the map. Inside the loop, the function motion_command is called to compute the pose of the robot.

Implement the missing parts in the file motion_command.m to compute the pose \mathbf{x}_t given \mathbf{x}_{t-1} and the odometry command \mathbf{u}_t . These vectors are in the following form:

$$\mathbf{x}_{t} = \begin{pmatrix} x \\ y \\ \theta \end{pmatrix}, \quad \mathbf{u}_{t} = \begin{pmatrix} \delta_{rot1} \\ \delta_{trans} \\ \delta_{rot2} \end{pmatrix},$$

where δ_{rot1} is the first rotation command, δ_{trans} is the translation command, and δ_{rot2} is the second rotation command. The pose is represented by the 3×1 vector x in motion_model.m. The odometry values can be accessed from the struct u using u.r1, u.t, and u.r2 respectively.

Compute the new robot pose according to the following motion model:

$$x_t = x_{t-1} + \delta_{trans} \cos(\theta_{t-1} + \delta_{rot1})$$

$$y_t = y_{t-1} + \delta_{trans} \sin(\theta_{t-1} + \delta_{rot1})$$

$$\theta_t = \theta_{t-1} + \delta_{rot1} + \delta_{rot2}$$

Test your implementation by running the main.m script. The script will generate a plot of the new robot pose at each time step and save an image of it in the plots directory. While debugging, run the program only for a few steps by replacing the for-loop in main.m by something along the lines of for t = 1:20. You can generate an animation from the saved images using avconv (sudo apt-get install libav-tools) or mencoder. With avconv you can use the following command to generate the animation from inside the plots directory:

avconv -r 10 -b 500000 -i odom_%03d.png odom.mp4

Exercise 3: Homogeneous transformations

A robot pose in a given frame is compactly represented as:

$$\mathbf{p} = \left(\begin{array}{c} x \\ y \\ \theta \end{array}\right)$$

Alternatively, homogeneous coordinates can be used. This simplifies transformations as they are modeled in matrix form:

$$\mathbf{T} = \begin{pmatrix} \mathbf{R} & \mathbf{t} \\ \mathbf{0} & 1 \end{pmatrix} = \begin{pmatrix} \cos \theta & -\sin \theta & x \\ \sin \theta & \cos \theta & y \\ 0 & 0 & 1 \end{pmatrix}$$

- (a) Implement two Octave funtions:
 - **v2t** takes as input the vector form of the robot pose and outputs the corresponding homogeneous transformation.
 - t2v takes as input a homogeneous transformation representing the robot pose in 2-D space and outputs the corresponding compact vector.

Test your implementation by chaining together four different transformations of your choice.

- (b) Given two robot poses $\mathbf{x}_1 = (x_1, y_1, \theta_1)^T$ and $\mathbf{x}_2 = (x_2, y_2, \theta_2)^T$, how do you get the relative transformation from \mathbf{x}_1 to \mathbf{x}_2 ?
- (c) A 2-D point (which we will often call "landmark" or "point of interest") in the plane can be expressed by its x and y values. In order to get a 3×1 vector and to multiply it with 3×3 matrices representing the robot poses, we need to express the landmark in homogeneous coordinates. This means adding a scaling value of 1:

$$\mathbf{poi} = \left(\begin{array}{c} x \\ y \\ 1 \end{array}\right).$$

Given a robot pose \mathbf{x}_t and an $\langle x, y \rangle$ observation \mathbf{z} of a landmark relative to \mathbf{x}_t ,

$$\mathbf{x}_t = \begin{pmatrix} 1 \\ 1 \\ \pi/2 \end{pmatrix}, \quad \mathbf{z} = \begin{pmatrix} 2 \\ 0 \end{pmatrix},$$

compute the location of the landmark.